SCHILLER

MEDICAL S.A.S.

FRED easy T

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1) GENERAL POINTS

- The TRAINING FRED device bears the CE CE-0459 mark in accordance with Directive 93/42/EEC relating to medical devices, and fulfils the fundamental requirements of Appendix I of that Directive.

- The device complies with standard EN 55011 class B.

- The safety of the patient and the user and trouble free device operation, shall only be guaranteed when original SCHILLER parts are correctly used. Further, only the accessories mentioned in these instructions for use or accessories that have been expressly authorised by SCHILLER may be used. If the device becomes damp as a result of temperature variations, it may only be put back into service once it is completely dry.

- The manufacturer shall not be held liable for the safety, reliability and performance of the device if:

- Assembly, extensions, adjustments, modifications or repairs are not performed by the manufacturer or by personnel authorised by the manufacturer.
- The electrical installation of the premises of use does not comply with locally applicable regulations.
- The device is not used in accordance with the instructions for use.

- The manufacturer may supply, on request, circuit diagrams, lists of components, descriptions, calibration instructions or any other information, for use by the qualified technical personnel of the user for repairing such parts of the device as have been stated to be repairable by the manufacturer. Such supply shall not in any case constitute permission or approval to modify or repair the device.

- This document describes a device type that complies with the safety standards applicable at the time of printing. SCHILLER reserves all rights to the circuits, processes, names, software and devices mentioned in these instructions for use.

- The user guide is an integral part of the device and should always be kept near the device. Close observance of the information given in the user guide is a prerequisite for using the device as intended for and correct operation and ensures patient and operator safety. Therefore, be sure to read the complete user guide.

- The user guide informs the device operator about the intended use, exact function, operation and required preventive maintenance. It is not a substitute for a product training.

- The quality assurance process used by SCHILLER is in compliance with international standards ISO 9001 and EN 46001.

- No part of these instructions for use may be reproduced without the written permission of SCHILLER.

- This device is not designed for any use other than that specifically described in this manual. Such use may be hazardous

MANUFACTURER:
SCHILLER Médical SAS
4, rue Louis Pasteur – BP 90050
F-67162 WISSEMBOURG Cedex

Tel: +33 (0)3.88.63.36.00
Fax: + 33 (0)3.88.94.12.82

E-mail: info@schiller.fr
Internet site: www.schiller-medical.com

Revision history:

Version	Date	Comment
01-01	24-04-2003	1st issue

2) INTENDED USE

The Training FRED is a simple cost-effective defibrillator simulator that is designed to be used exclusively for training. In no event may it be used as a conventional defibrillator. It has been designed for the purpose of training users of Schiller's FRED Easy semi-automatic defibrillator. The device provides several simulations or scenarios (from 0 to 8) that aim to help beginners to familiarise themselves with the FRED Easy and show that they have the knowledge required to use it in case of emergency. The device also allows the trainer to work using an infrared remote control so that he can intervene in real time on the scenario in progress.

3) CHARACTERISTICS

The FRED easy Trainer allows the **simulated** delivery of electric shocks. There is no high voltage in the device, and training is risk free.



The FRED Easy trainer may be used on all types of mannequins as it does not deliver an electric shock. It may only be used with the training electrodes sold by Schiller.


The FRED Easy Trainer may be used with an infrared remote control. The remote control allows the trainer to control the training actions.

The FRED Easy Trainer has 9 pre-recorded training scenarios which simulate realistic cardiac arrest configurations.

The basic equipment of the Training FRED consists of:

- 4 LR6 batteries which may be replaced by rechargeable batteries.
- training electrodes consisting of a single use electrode part and a reusable cable. The cable is fitted with a magnet allowing the device to recognise the presence of the electrodes in the device without any electrical contact.

	Reusable training electrode with adhesive disk for fixing the reusable cable.
	Reusable cable fitted with a magnet in the connector so that the FRED Easy T knows the electrodes are present.

The Training FRED offers the user a choice of 9 possible scenarios numbered 0 to 8, which are described in the paragraph below.

4) SCENARIO SELECTION PROCEDURE


List of scenarios:

SNR= Shock not recommended

SR = Shock recommended

Elect prob = Electrode problem

Scenario number 0	SNR
Scenario number 1	SR/SR/SNR
Scenario number 2	SR/SR/SR/SNR
Scenario number 3	SR/SR/SR/SR/SNR
Scenario number 4	SR/SR/Elect prob/SR/SNR
Scenario number 5	Elect prob/SR/SR/SNR
Scenario number 6	SR/SNR/SR/SNR
Scenario number 7	SNR/SR/SNR/SNR
Scenario number 8	SR/SR/SR
Scenario number 9	Infrared control function

1) The device is initially completely off, the green diode (3) is not flashing.

2) Press the green key (1) to switch on the device, the device comes on on the last selected scenario.

3) Place the device on standby mode by pressing the green key (1) for a few seconds.

4) At the same time press the orange key (3) first and the green key (1) to access the configuration menu.

5) Select automatic or manual analysis mode according to the configuration of the FRED Easy which you wish to use for training. Select by pressing successively on the green key (1) and change to the next menu by pressing on the orange key (3)

6) In the same way, select the PAD configuration (**P**ublic **A**ccess **D**efibrillator):

none = Configuration without information messages (none of the messages 1 to 7 in the table in Appendix A will be displayed).

Short = Configuration with brief information messages (the messages 1 and 7 in the table in Appendix A will be displayed)

True = Configuration with full information messages (all the messages 1 to 7 in the table in Appendix A will be displayed)

7) In the same way, select your scenario 0 to 8 or 9 for infrared control mode.

8) The next menu proposes that you save your selection. Validate by pressing the green key (1). The device then switches off and is placed on standby mode. The green diode (3) flashes to inform the user that the device is operating normally.

The user is then ready to start the device normally and it will operate according to the scenario selected.

In the case of the selection of scenario number 9 the device awaits the trainer's orders before changing status. These orders are transmitted using an infrared control and received in the device by means of a receiver situated as shown (4) in the photograph above.

The functions associated with keys on the infrared control are as follows:

Key 0	Immediate move to the next step
Key 1	‘Support with dealing with the victim’ information messages
Key 2	Shock recommended
Key 3	Shock not recommended
Key 4	The rhythm has changed
Key 5	Press the electrodes (Contact problem)
Key 6	Switch off the green diode which means device is off
Key 7	Switch off flashing green diode
Key 8	Check the circulation
Key 9	Pause


The device is only totally switched off when the green diode is no longer flashing. Do not leave it on standby, as this wears the batteries out unnecessarily.

To switch off the device:

Whilst a scenario is in operation, press at the same time for a few seconds on the orange key and the green key until the device switches off completely. It is no longer consuming energy.

Example of a scenario:

Scenario number 2

5) SAFETY INSTRUCTIONS

The following points shall be complied with in particular:

- i) The device has been designed for training purposes **only**. It may not in any event be used for defibrillation.
- j) The device is used with training electrodes.
- k) These instructions for use describe the complete FRED Easy training device with all its optional features.
- l) The device is not designed for any use other than that described in these instructions. Such use is forbidden.
- m) The user shall be entirely liable for use of accessories of makes other than SCHILLER.
- n) The user must be familiar with the functioning of the device.
- o) Before each use, the user must ensure the device is functioning correctly.
- p) The device is not designed for use in areas where an explosion hazard may occurs

6) CLEANING AND MAINTENANCE

The outer surfaces of the device may be cleaned with cotton wool moistened with alcohol.

Important

**Risk of damage to the device –
Never use phenol based cleaning agents or agents containing peroxide
derivatives to disinfect the surfaces of the device housing**

During cleaning, the recyclable batteries must be removed from the device. No liquid may be allowed to enter the device. If liquid does penetrate into the device, clean the device and have it verified completely.

7) SERVICING AND DISPOSAL AT THE END OF ITS SERVICE LIFE

Have device performance checked once a year. At the end of its service life, the device and the accessories must be disposed of in compliance with the local regulations. Apart from the plug-in batteries, the device does not contain hazardous material and can be disposed of like any other electronic equipment.

8) SPECIFICATIONS

- Operating mode:

Training in semi-automatic defibrillation

- Power supply:

4 LR6 batteries

- 1 to 100% on-screen battery indicator corresponding to the battery voltage range from 4.5 volts minimum to 6 volts.

- Electrode and battery fault simulation
- Detection of electrode presence by means of a magnet placed in the connector of the training electrodes.
- Connector for the training electrodes


PHYSICAL AND ENVIRONMENT REQUIREMENTS

- Transport /storage:
temperature $-30...+50^{\circ}\text{C}$
Relative humidity rate 30...95 %, without condensation
Atmospheric pressure 500...1060hPa
- Service:
Temperature $0...+50^{\circ}\text{C}$
Relative humidity rate 30...95 %, without condensation
Atmospheric pressure 700...1060hPa

DIMENSIONS AND WEIGHT

Width	230 mm
Depth	220 mm
Height	70 mm
Weight	<1 kg

9) DISPLAY AND CONTROLS

10) WARRANTY

Schiller Médical guarantees its devices to the purchaser against all material or manufacturing defects for a period of one year as from the original purchase date. During this one-year period, Schiller Médical, as soon as it receives a device presenting a material or manufacturing defect returned accompanied by written notification of the defect, will repair or replace the defective parts or the complete product if necessary.

All postal, shipping or handling costs will be borne exclusively by the user.

The manufacturer shall supply, on request, circuit diagrams, lists of components, descriptions, calibration instructions or any other information, for use by the qualified technical personnel of the user for repairing such parts of the device as have been stated to be repairable by the manufacturer. Such supply shall not in any case constitute permission or approval to modify or repair the devices.

All rights are reserved for the circuits, processes and names mentioned in these instructions for use.

Schiller Médical shall not be held liable as a result of this warranty for any secondary or indirect damage, or in the event of unauthorised repairs or modifications having been made or undertaken, or when the whole or a part of the product has been damaged by accident, incorrect use or deliberate misuse. This warranty does not cover the batteries, normal wear and tear, marks or any other irregularity of appearance which does not prevent or hinder the performance of the device.

11) Information Order

Ref	Designation
53108	4 Batteries LR6
	Path electrodes for training

APPENDIX A

1	Spoken text	Text displayed
	<u>Basic Life Support for full version</u>	
1.	Stay calm, Fred is here to help your	
2	Shake the patient and ask him/her "How are you ? What happened ?"	<u>SHAKE THE PATIENT</u>
3	If there is no reaction, seek help and have someone call the emergency services	<u>CALL THE EMERGENCY SERVICES</u>
4	Check if the patient is breathing and if he/she is showing any signs of circulation	<u>CHECK BREATHING AND FOR SIGNS OF CIRCULATION</u>
5	If the patient is breathing, roll him/her over onto his/her side	<u>BREATHING – PATIENT ONTO SIDE</u>
6	If the patient is not breathing, but is showing signs of circulation, administer artificial respiration until the emergency services arrive	<u>NO BREATHING – ARTIFICIAL RESPIRATION</u>
7	If the patient is not showing any signs of circulation, stick the electrodes on his/her bare chest as indicated	<u>NO SIGNS OF CIRCULATION – STICK ON THE ELECTRODES</u>

Step 4 – 7 repeat (every 20 sec) until FRED recognises that the electrodes are connected to the patient. Then continue with steps 8, or 9 depending on the analysis mode (manual or automatic) or with step 12 in the event of asystole.

8	For analysis, press the green key	<u>FOR ANALYSIS, PRESS THE GREEN KEY</u>
9	Do not touch the patient, analysis in progress	<u>DO NOT TOUCH THE PATIENT- ANALYSIS IN PROGRESS</u>

SHOCK RECOMMENDED

10	Shock recommended. Charge. Do not touch the patient	<u>CHARGE-DO NOT TOUCH THE PATIENT</u>
11	Administer the shock. Do not touch the patient. Press the orange shock key	<u>SHOCK: PRESS THE ORANGE SHOCK KEY</u>

Automatic repetition of steps 9 to 11 twice – if fibrillation persists, then

12	Check if the patient is breathing and if he/she is showing signs of circulation	<u>CHECK BREATHING AND SIGNS OF CIRCULATION</u>
13	No signs of circulation: alternately administer artificial respiration twice and 15 cardiac massages. If he/she is showing signs of circulation, roll the patient onto his/her side	<u>IF NONE : 2 X ARTIFICIAL RESPIRATION & 15 MESSAGES IF SIGNS OF CIRCULATION - PATIENT ONTO SIDE</u>

One minute after step 13? the device goes back to step 8 or 9 depending on the analysis mode (manual or automatic)

OR: SHOCK NOT RECOMMENDED

14	Shock not recommended	
----	-----------------------	--

The device goes to step 12.

IF NORMAL RHYTHM AFTER SHOCK RECOMMENDED

15		THE RHYTHM HAS CHANGED
----	--	-------------------------------

The device goes to step 12.

Technical fault:

16	Press the electrodes	PRESS THE ELECTRODES
17	Stick the electrodes on the patient's bare chest as indicated	STICK ON THE ELECTRODES
18	Movements detected. Stop movements	MOVEMENTS DETECTED. STOP MOVEMENTS